

## Accessories

Photo	Name	Ordering Code	Size: ( L x W ) mm
	8-chamber Arm Sleeve	220081	813.41 x 302.14
	8-chamber Leg Sleeve	220082	S: 866.15 x 373.95
		220082M	M: 905.78 x 444.14
		220082L	L: 948.68 x 490.25
	4-chamber Waist Sleeve	220083	M: 1360.00 x 400.00
		220083L	L: 1600.00 x 400.00
	8-chamber Full Body Cuff	220016FB	1463.43 x 611.14

## Certification



National Medical Products Administration



CE(MDD)



ISO13485 / MDSAP

# Longest

Guangzhou Longest Science & Technology Co., Ltd.

Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, China 511399

Tel: 86-20-66353999 Fax: 86-20-66353920

E-mail: export@longest.cn Website: www.gzlongest.com



# LGT-2200WM

## Compression Therapy Device

Health Life Longest Care

# Longest






# LGT-2200WM


## Compression Therapy Device

GT-2200WM is a sequential pneumatic pressure system that through program controlled inflation and deflation of the sleeves covered in arms, legs, waist, and back to realize therapeutic effect



### ⊕ Specifications

- ⊗ 8 chambers
- ⊗ 8 compression modes
- ⊗ Time range: 1-99 mins
- ⊗ Pressure range: 40-190 mmHg
- ⊗ Interface: Touch screen
- ⊗ Weight: 9 kg
- ⊗ Size: 366mm(W) × 393mm(L) × 194mm(H)



### ⊕ Mechanism and functions of LGT-2200WM

#### ⊗ For physio-

The inflation and deflation of the sleeves enhances the backflow of venous blood and lymph, which can improve the blood circulation of diseased areas. This will lead a series of beneficial effect such as improving the function of peripheral blood vessels, removing the edema, accelerating the healing process, and preventing Muscular Dystrophy.

#### ⊗ For Aesthetics-

The Multi Lymphatic drainage system, LGT-2200WM is used for the purpose of drainage of the dissolved or melted fat cells generated after various obesity equipments by improving lymph glands or blood vessels that constricted or damaged due to cellulite or excessive fat cells on the lymphatic drainage step, one of general obesity lipolysis steps.



### ⊕ Features

- ⊗ User-friendly interface, just a few steps to finish your setting
- ⊗ Gentle inflation and deflation
- ⊗ Real-time pressure, time, and inflation part monitoring (visible in the interface)
- ⊗ Durable & eas - use sleeves
- ⊗ Pressure unit switchable between mmHg and kPa
- ⊗ 2 pcs full body cuffs can be connected with one device at the same time

### ⊕ Application

- ⊗ Lymph drainage
- ⊗ Blood circulation
- ⊗ Detoxification
- ⊗ Cellulite treatment
- ⊗ Body shaping in legs, abdomen, buttocks, and arms
- ⊗ Strengthening of loose skin
- ⊗ Relaxation of stress
- ⊗ Relieve of fatigue